



AN 402

## Strained Si - Surface Morphology Measured by AFM

May 7, 2007 (Version 3.0)

### Discussion


Enhanced channel mobility of CMOS devices has become a key technology goal to continue Moore's Law. One promising technical approach is to use a strained Si channel created by epitaxially depositing a thin (~ 10 nm) silicon layer on top of a strain-relaxed  $\text{Si}_{1-x}\text{Ge}_x$  buffer layer film which has been grown on top of a silicon substrate. The difference in lattice constant between the two films results in a tensile strain in the top thin silicon layer.

However, the growth of the strain-relaxed  $\text{Si}_{1-x}\text{Ge}_x$  layer on top of the silicon substrate also results in the creation of many dislocations which then create defects in the thin silicon channel layer. Usually, two types of dislocations, misfit and threading, are observed in the strain-relaxed  $\text{Si}_{1-x}\text{Ge}_x$  layer.


Because the presence of threading dislocations can seriously degrade the device performance, graded  $\text{Si}_{1-x}\text{Ge}_x$  layers are grown to minimize the density of threading dislocations. Furthermore, the misfit dislocations will lead to crosshatch patterns on the top surface of  $\text{Si}_{1-x}\text{Ge}_x$  layer. When the thin tensile-strained Si layer is epitaxially grown on this strain-relaxed  $\text{Si}_{1-x}\text{Ge}_x$  surface, the cross-hatch patterns are conformal.

Surface height variations of the crosshatch patterns can be significantly large as compared to the thickness of the tensile-strained Si layer. It is therefore important to minimize the surface variation (crosshatch patterns) of the tensile-strained Si surface.

For this purpose, AFM can provide quantitative surface roughness measurements and mapping of the crosshatch dimensions to assist in the process development and monitoring. Figures 1 and 2 show AFM images of two different tensile-strained Si surfaces. Sample 2 is smoother than sample 1 (0.6 nm compared to 1.7 nm), and the crosshatch dimensions (spatial frequencies) are very different (Neither sample received a CMP step to smooth the surface.)



**Figure 1. Tensile-strained Si surface of sample 1 Rrms: 1.7nm**



**Figure 2. Tensile-strained Si surface of sample 2 Rrms: 0.6nm**

---

**United States Locations**

Tempe, Arizona  
+1 480 239 0602 info.az@eaglabs.com  
+1 602 470 2655 fax

Sunnyvale, California  
810 Kifer Road  
+1 408 530 3500 info.ca@eaglabs.com  
+1 408 530 3501 fax

1135 E Arques Avenue  
+1 408 738 3033  
+1 408 738 3035 fax

785 Lucerne Drive  
+1 408 737 3892  
+1 408 737 3916 fax

Peabody, Massachusetts  
+1 978 278 9500 info.ma@eaglabs.com  
+1 978 278 9501 fax

Chanhassen, Minnesota  
+1 952 828 6411 info.mn@eaglabs.com  
+1 952 828 6449 fax

East Windsor, New Jersey  
+1 609 371 4800 info.nj@eaglabs.com  
+1 609 371 5666 fax

Syracuse, New York  
+1 315 431 9900 info.ny@eaglabs.com  
+1 315 431 9800 fax

Raleigh, North Carolina  
+1 919 829 7041 info.nc@eaglabs.com  
+1 919 829 5518 fax

Round Rock, Texas  
+1 512 671 9500 info.tx@eaglabs.com  
+1 512 671 9501 fax

**International Locations**

Shanghai, China  
+ 86 21 6879 6088 info.cn@eaglabs.com  
+ 86 21 6879 9086 fax

Tournefeuille, France  
+ 33 5 61 73 15 29 info.fr@eaglabs.com  
+ 33 5 61 73 15 67 fax

Frankfurt, Germany  
+ 49 (0) 693053213 info.de@eaglabs.com  
+ 49 (0) 69307941 fax

Tokyo, Japan  
+ 81 3 5396 0531 info.jp@eaglabs.com  
+ 81 3 5396 1930 fax

HsinChu, Taiwan  
+ 886 3 5632303 info.tw@eaglabs.com  
+ 886 3 5632306 fax

Uxbridge, United Kingdom  
+ 44 (0) 1895 811194 info.uk@eaglabs.com  
+ 44 (0) 1895 810350 fax