

Analytical Solutions

May 9, 2007 (Version 2.0)

for BioTechnology


BN 1423


XPS Analysis of Disposable Gloves

Discussion

Many types of processing involve some degree of sample handling. In industries where cleanliness is critical to quality, disposable elastic gloves are widely used to protect critical surfaces from inadvertent contamination. If you are relying on your gloves to provide a clean barrier between skin and a critical surface, you must be confident that the barrier itself is not a problem. X-ray photoelectron spectroscopy (XPS), also called Electron Spectroscopy for Chemical Analysis (ESCA), provides a sensitive tool for quantitatively evaluating the cleanliness of the gloves that you depend on. Below are results from the XPS analysis of five different brands of disposable elastic gloves. The table shows atomic percentages of all the elements detected on the five gloves. It reveals that varying levels of Si (as silicone), S, Cl, Ca and Zn are detected on the gloves, in addition to the expected polymer components.

Atomic % of Elements Detected on Gloves								
	C	N	O	Si	S	Cl	Ca	Zn
Glove 1	85	1.5	8.6	1.5	0.5	2.6	0.2	0.5
Glove 2	78	3.3	10	-	1.1	6.0	-	1.2
Glove 3	54	0.3	26	19	-	0.1	0.5	0.1
Glove 4	76	2.3	13.2	-	1.4	4.9	0.5	1.7
Glove 5	71	0.1	17	12	-	0.1	-	0.1





United States Locations

Tempe, Arizona
+1 480 239 0602 info.az@eaglabs.com
+1 602 470 2655 fax

Sunnyvale, California
810 Kifer Road
+1 408 530 3500 info.ca@eaglabs.com
+1 408 530 3501 fax

1135 E Arques Avenue
+1 408 738 3033
+1 408 738 3035 fax

785 Lucerne Drive
+1 408 737 3892
+1 408 737 3916 fax

Peabody, Massachusetts
+1 978 278 9500 info.ma@eaglabs.com
+1 978 278 9501 fax

Chanhassen, Minnesota
+1 952 828 6411 info.mn@eaglabs.com
+1 952 828 6449 fax

East Windsor, New Jersey
+1 609 371 4800 info.nj@eaglabs.com
+1 609 371 5666 fax

Syracuse, New York
+1 315 431 9900 info.ny@eaglabs.com
+1 315 431 9800 fax

Raleigh, North Carolina
+1 919 829 7041 info.nc@eaglabs.com
+1 919 829 5518 fax

Round Rock, Texas
+1 512 671 9500 info.tx@eaglabs.com
+1 512 671 9501 fax

International Locations

Shanghai, China
+ 86 21 6879 6088 info.cn@eaglabs.com
+ 86 21 6879 9086 fax

Tournefeuille, France
+ 33 5 61 73 15 29 info.fr@eaglabs.com
+ 33 5 61 73 15 67 fax

Frankfurt, Germany
+ 49 (0) 693053213 info.de@eaglabs.com
+ 49 (0) 69307941 fax

Tokyo, Japan
+ 81 3 5396 0531 info.jp@eaglabs.com
+ 81 3 5396 1930 fax

HsinChu, Taiwan
+ 886 3 5632303 info.tw@eaglabs.com
+ 886 3 5632306 fax

Uxbridge, United Kingdom
+ 44 (0) 1895 811194 info.uk@eaglabs.com
+ 44 (0) 1895 810350 fax