



Dose and Oxide Thickness Determination by XPS of Plasma-Implanted P in Si

Wesley Nieveen
Evans Analytical Group

Poster presented at USJ 2005
Daytona Beach, FL
June 5-8, 2005

United States Locations

Tempe, Arizona

+1 480 239 0602
info.az@eaglabs.com

Sunnyvale, California

+1 408 530 3500
+1 408 738 3033
+1 408 737 3892
info.ca@eaglabs.com

Peabody, Massachusetts

+1 978 278 9500
info.ma@eaglabs.com

Chanhassen, Minnesota

+1 952 828 6411
info.mn@eaglabs.com

East Windsor, New Jersey

+1 609 371 4800
info.nj@eaglabs.com

Syracuse, New York

+1 315 431 9900
info.ny@eaglabs.com

Raleigh, North Carolina

+1 919 829 7041
info.nc@eaglabs.com

Round Rock, Texas

+1 512 671 9500
info.tx@eaglabs.com

International Locations

Shanghai, China

+ 86 21 6879 6088
info.cn@eaglabs.com

Tournefeuille, France

+ 33 5 61 73 15 29
info.fr@eaglabs.com

Frankfurt, Germany

+ 49 (0) 693053213
info.de@eaglabs.com

Tokyo, Japan

+ 81 3 5396 0531
info.jp@eaglabs.com

HsinChu, Taiwan

+ 886 3 5632303
info.tw@eaglabs.com

Uxbridge, United Kingdom

+ 44 (0) 1895 811194
info.uk@eaglabs.com

Specialists in Materials Characterization

INTRODUCTION

Future Node Implant Issues:

- What are the difficulties in measuring dopants?
 - ITRS roadmap has junction depth (X_j) of <10nm at 45nm node
 - Lateral abruptness, channel length, resistivity, and other ITRS targets are pushing dopant levels to $>E+20/cm^3$
 - Dopant doses required are $>$ mid-E+14 atoms/cm².
- Two often cited methods for ULE/USJ implant to provide high dose at shallow controlled depths: [1]
 - (Low Energy) plasma implant using gas phase species, e.g. PH₃, BF₃, B₂H₆ etc.
 - 'Standard' implantation at low energy with poly-atomic species, e.g. B₁₀H₁₈, B₁₈H₂₂, P_xH_{3x} etc.
- Activation of implant still required – methods under investigation include:
 - 'spike' anneals (RTA, laser, etc.)
 - solid phase epitaxial re-growth


¹W. J. Taylor et al., *Materials Challenges for CMOS Junctions*, Mat. Res. Soc. Symp. Proc. Vol. 810, 2004

Some Fundamental Problems Measuring ULE/USJ Implant Dopants

- At $X_j < 15nm$ and Dose $> E+20$ atoms/cm³
 - ⇒ Concentration is on the order of a few at%
 - ∴ This is ideally suited to measurement by XPS
- Regardless of technique used to implant dopant, some amount of surface oxide will be present on the Si.
 - X% of dopant dose will be implanted into the Si substrate.
 - Y% of dopant dose will be implanted in SiO₂ layer.
 - Z% of dopant dose may be lost or unaccounted for in implant process.
- Anneal/Activation may re-distribute dose in oxide, dose in substrate, and/or cause additional loss of implant dose.

BACKGROUND


Focus on the Oxide Layer Problem



Junction Depth, X_j :
If we define X_j as the depth at which the implant attains a concentration of 1×10^{18} atoms/cm³, **Then**, inclusion or exclusion of the oxide layer can change X_j by $>70\%$ in this case!

This also implies that the dose in the oxide must be known separately and independently from the dose in the Si substrate.


XPS Perspective of Plasma Doped P in Si



No PO_x atoms are in substrate and no P⁰ atoms are in oxide.

XPS Assumptions & Caveats

- XPS assumes: All the dopant is within the depth of information (DOI). For P in SiO₂/Si, all P must be within ~150Å of surface (oxide layer + Si substrate)
- XPS assumes: All the dopant is uniformly distributed within the oxide layer and the substrate layer.
- If P distribution within a layer (oxide or Si) is peaked toward surface, XPS will overestimate P dose.
- If P distribution within a layer is peaked toward substrate, XPS will underestimate P dose.



SAMPLES

Plasma Doped P in Si Samples

- Samples supplied by a major international research institute
- Plasma method of implantation of P into Si (PH₃ in carrier gas)
- Samples consist of:
 - As-implanted wafers (no post plasma treatment)
 - Annealed wafers at 900C for 10 seconds.
- Portions of 200mm wafers cleaved and analyzed by XPS and SIMS
- P present on as-implanted and annealed samples. Other 'impurities' (probably from plasma implanter) also present.

RESULTS & DISCUSSION

Measured XPS Oxide Thickness (T_{ox}) and Calculated P Dose in Oxide (P_{ox}) and P Dose in Silicon Substrate (P_{sub})

	T_{ox} (Å)	P_{ox} Dose (at/cm ²)	P_{sub} Dose (at/cm ²)	Total P Dose (at/cm ²)
P in Si; as-implanted	43.87	4.71E+15	1.92E+15	6.63E+15
P in Si; RTA @900C 10sec	42.96	1.07E+15	1.42E+15	2.48E+15

Intermediate Conclusion:


- The following observations can be made:
 - As-implanted: P dose in oxide $>$ P dose in substrate.
 - After RTA 900C-10s: P dose in oxide $<$ P dose in substrate.
 - P dose as-implanted in oxide $>$ P dose after RTA anneal in oxide
 - P dose as-implanted in substrate \approx P dose after RTA anneal in substrate

What happened to the P?

- As-implanted: most of P in oxide w/ small amount in substrate.
- RTA anneal @ 900C for 10s: more P in substrate after anneal.

	Oat%	Siat%	Total P
As-implanted	67.4	24.1	6.63E+15/cm ²
RTA @900C 10sec	64.9	31.8	2.48E+15/cm ²

Where did the rest of the P (~4.2E+15 atoms/cm²) go after anneal?



'As-implanted' vs. RTA anneal @ 900C-10s change to oxide layer.

Thickness of oxide ~ identical before & after anneal.

	T_{ox} (Å)	P_{ox} Dose (at/cm ²)	P_{sub} Dose (at/cm ²)
As-implanted	43.87	4.71E+15	1.92E+15
RTA @900C/10s	42.96	1.07E+15	1.42E+15
Delta	-0.91 (-2.07%)	-3.64E+15 (-77.28%)	-0.50E+15 (-26.04%)


Intermediate Conclusion:

- If oxide thickness did not change from as-implanted to annealed; **Then** ~60% loss of P is due to:
 - P evaporating during anneal OR
 - P migrating to depth $>$ DOI (XPS)

(DOI = Depth of Information \approx 160Å @ 75° TOA)

Is XPS measuring all the dopant? From SIMS Profiles:

- XPS 'sees' $>99.9\%$ of P! (All of the P is within the XPS depth of information.)
- P concentration peaks just below the surface and then decays steadily.
- Appears highest concentration of P is at surface with deeper sub-surface P peak compared to as-implanted sample.



RESULTS & DISCUSSION


SIMS-XPS Results P Plasma Implant in Si:

- SIMS shows all the P is detectable by XPS within top 100-150Å!
- Results imply some P is lost to vacuum/air during RTA.
- Shows some P is 'driven' into substrate after anneal.
- XPS cannot accurately determine X_j – requires SIMS profile. *Two-technique measurement gives accurate and precise USJ/ULE dose and X_j values.*

XPS Precision and Accuracy for ULE P Dose Measurements in Si

- XPS repeatability (precision) is very high (based on SiON experience)
 - Requirements:
 - Sample stability
 - Stable achromatic/monochromatic source
 - No 'flood' gun/neutralization
 - High signal-to-noise
 - Short measurement time
 - Minimal or no operator intervention
 - Robust algorithm
 - High stability Instrument –stage, electronics, analyzer, vacuum, etc.
 - Long term reproducibility
 - 'Clean' samples
 - Accurate & calibrated – tool matching
 - Achievements of SiON HPM (High Precision Measurement)
 - $4\text{Å} \geq T_{ox} \geq 150\text{Å}$
 - ($\geq 90\text{Å}$ with reduced precision)
 - T_{ox} 'static' precision (1σ) $< 0.05\text{Å}$
 - T_{ox} typical precision $< 0.2\text{Å}$
 - N Dose (min.) $> 1 \times 10^{14}$ at/cm² (lower w/ decreased precision)
 - N Dose: Rel. Std. Dev. $< 1.0\%$, typical precision $< 0.05 \times 10^{15}/cm^2$
 - XPS HPM (High Precision Measurements) are highly specific and require specific instrument configurations and tool matching.

- Absolute accuracy will depend somewhat on P distribution within layer.
- Uniformity of:
 - Dose in substrate
 - Dose in oxide
 - Oxide thickness
 across a wafer will be critical.
- Dose (in oxide or substrate) does not have to correlate to oxide thickness: therefore it is important to measure both independently



CONCLUSIONS

- A single XPS measurement can provide:
 - accurate and precise oxide layer thickness.
 - total P dose and P dose of oxide and substrate regions individually.
- Excellent short term and long term precision of XPS measurement (demonstrated by SiON measurements).
- XPS should have good sensitivity to process changes and conditions.
- XPS can be used for P dose (radial) implant uniformity.

Appendix:

Uniform Constant Concentration Assumption: % Error for Sloped Line +/- Distribution

